Istruzioni per l'uso
Instructions for use
Instructions d'utilisation
Bedienungsanleitung
Οδηγίες χρήσης

TERMOFORO

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
DESCRIPTION OF THE APPLIANCE

- | | |
|--|--|
| 1. Heated area | 5. Opening of integrated pocket for power bank |
| 2. ON/OFF and temperature selection button | 6. Product label |
| 3. Temperature visualisation LED | 7. USB-A connection cable |
| 4. Integrated pocket for power bank | 8. Keyhole for closing the scarf |

TECHNICAL DATA

TYPE S6101 5 V

TENACTA GROUP S.p.A.
VIA PIEMONTE 5, IT
24052 AZZANO S.P.

Dear customer, thank you for choosing our product. We feel certain you will appreciate its quality and reliability as it was designed and manufactured with customer satisfaction in mind. These instructions for use are compliant with the European Standard EN 82079.

INDEX

Illustrative guide	page 15
Description of the appliance	page 15
Technical data	page 15
Symbols	page 16
Important instructions	page 17
Use	page 21
Cleaning and maintenance	page 23
Storing the heating scarf	page 24
Troubleshooting	page 25
Disposal	page 26
Support and warranty	page 26

SYMBOLS

Direct current	Prohibited
Carefully read the instruction manual	Important note
Warning	Not suitable for children under the age of 3
Do not insert pins	Do not bleach
Wash at 30°C	Do not dry-clean
Do not iron	Do not tumble dry

IMPORTANT INSTRUCTIONS. RETAIN FOR FUTURE USE.

WARNING! Instructions and warnings for safe use.

Before using this appliance, carefully read the instructions for use and in particular the safety notes and warnings, which must be complied with. Keep this manual together with its illustrative guide for future consultation. Should you transfer the appliance to another user, make sure you hand over this documentation as well.

NOTE: should any part of this manual be difficult to understand or any doubt arise, contact the company at the address indicated on the last page before using the appliance.

NOTE: persons with impaired vision may consult these instructions for use in digital format at the website www.tenactagroup.com

- After unpacking, make sure the appliance is intact, complete with all its parts as shown in the illustrative guide and with no sign of damage. If in any doubt, do not use the appliance and contact an authorised service centre.

WARNING! Risk of suffocation.

Children shall not play with the packaging. Keep the plastic bag out of the reach of children.

- This appliance must be used only for the purpose for which it was designed, namely as a pad for home use. Any other use is considered not compliant with the intended use and therefore dangerous.

Children under the age of three are not to use this appliance due to their inability to react to overheating.

- This appliance can be used by children older than 3 years and younger than 8 years under supervision and with the control always set to minimum temperature (temperature turquoise).

- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.
- This appliance must not be used by persons insensitive to heat and other very vulnerable persons who are unable to react to overheating.
- Always disconnect the USB connection cable from the power bank before carrying out any cleaning or maintenance, and also when the appliance is not being used.
- Examine the appliance frequently for signs of wear or damage. If there are such signs, if the appliance has been misused or does not work, contact an authorised technical service centre before switching it on again.
- If the supply cord is damaged, take the whole appliance to an authorised technical service centre so as to prevent all risks.
- This appliance is not intended for medical use in hospitals.

-  DO NOT use if wet.

DO NOT INSERT PINS OR NEEDLES.

 When operating the appliance over a prolonged period of time **DO NOT** exceed the temperature levels stated for continuous usage.

If the appliance is set and used at a temperature higher than those indicated for continuous usage over an extended period of time, the user may suffer from skin burns or heat-related injuries.

- When storing the appliance, allow it to cool down before folding.

WARNING! DO NOT crease the appliance by placing items on top of it during storage.

WARNING! Follow indications in the relevant chapters of this manual on how to use, clean, maintain or store the appliance.

- **The use of creams or ointments combined with the heat of the appliance may increase the risk of skin burns.**
- **Persons with a pacemaker must NOT use this appliance.**
- **⊘ DO NOT use the appliance with wet hands, damp feet or barefoot.**
- **⊘ DO NOT pull or lift the appliance by the supply cord.**
- **If the appliance is faulty or malfunctions, switch it off and do not tamper with it. Contact an authorised technical service centre for any repairs.**
- **This appliance shall not come into contact with liquids, sprays, steam. For any cleaning and maintenance operation refer to the instructions in the relevant chapter.**
- **⊘ DO NOT use this appliance to warm animals.**
- **See the external packaging for the characteristics of the appliance.**

USE

- The product starts working when connected to a power bank.
- The characteristics of the power bank must meet at least the following technical requisites.
5000mAh capacity
1 USB-C recharging port Input DC5V, 2A
1 micro-USB recharging port Input DC5V, 2A IT
1 USB-A recharging port Output DC5V, 2A

- The USB-A connection cable (7) is compatible with a USB-A to USB-C adapter

CONNECTING THE POWER BANK TO THE HEATED SCARF, AND SWITCHING ON

- Connect the USB connection cable supplied (7) to the USB-A connection port on the power bank. The temperature visualisation LED (3) will be red
- Insert the power bank in the integrated pocket (4) through the relative opening (5)
- Open out the scarf completely and then wrap it around your neck.
- Wear the scarf with the fleece on the inside, next to your neck.
- Thread the right end of the scarf through the keyhole (8) in the left end. (FIGURE A)
- Close and tighten as preferred. (FIGURE B)

CHOOSING THE TEMPERATURE

- Press the ON/OFF and temperature selection button (2) to set the required temperature. There are 3 settings:
OFF: temperature visualisation LED (3) blue
Low: temperature visualisation LED (3) turquoise
Medium: temperature visualisation LED (3) green
High: temperature visualisation LED (3) red
- To heat the scarf quickly, use the maximum temperature - «red». You can then set a lower temperature level, depending on the degree of heat required.
- Pressing the ON/OFF and temperature selection button (2), the temperature changes in the sequence High Red → Medium Green → Low Turquoise → High Red → Medium Green → etc.
- If the scarf is used for a long time, set the low temperature level - «turquoise».

AUTOMATIC SWITCH-OFF FUNCTION

- The scarf has an automatic switch-off system. After it has been working for about 90 minutes, the temperature visualisation LED (3) turns blue to indicate that the scarf has switched off automatically.
- To put the scarf on standby, press and hold the ON/OFF and

temperature selection button (2) for about 3 seconds. The temperature visualisation LED (3) will be blue.

- To reactivate the scarf, press the ON/OFF and temperature selection button (2).

SWITCHING OFF

- When you have finished using the scarf, disconnect the USB connection cable (6) from the USB-A connection port of the power bank.

CAUTION! Disconnect the power bank when the scarf is not being used.

- If the power bank discharges completely during use, the scarf heating function will automatically stop working.

CLEANING AND MAINTENANCE

CAUTION! Hang the appliance without applying pegs so as not to damage it.

WARNING! Always remove the power bank before carrying out any cleaning or maintenance.

WARNING! DO NOT use the device if damp

BEFORE WASHING

- Disconnect the power bank
- Disconnect the USB connection cable (7) from the power bank connection port.

CLEANING THE SCARF

- The scarf can be washed in lukewarm water (30°C). The use of a delicate detergent is recommended.

CLEANING THE SCARF IN THE WASHING MACHINE

- The scarf can machine washed at 30°C.
- Use only the delicate cycle for spinning. DO NOT run an additional spin cycle.
-  DO NOT iron the scarf.

DRYING IN THE OPEN AIR

- Spread the scarf out horizontally on a clothes dryer and leave to dry completely.

WARNING! Dry the appliance flat without using pegs, to avoid damaging it. Let it dry well, laid out flat on a clothes dryer.

-  DO NOT dry clean.
- DO NOT wring it out.
- DO NOT tumble dry.

STORING THE HEATED SCARF

When not in use, store as follows:

- If the scarf has just been washed, wait until it is fully dry.
- For better storage, it is advisable to place the scarf in the packaging provided at the time of purchase.
- When putting the scarf away, allow it to cool down before folding it.
- Keep the scarf in a dry place.

Examine the scarf frequently for signs of wear or damage before using it. If you notice such signs, or any malfunctioning, contact a service centre.

TROUBLESHOOTING

This chapter looks at the most common problems relating to the use of the appliance. If the problems cannot be solved with the information below, please contact the Authorised Assistance Centre.

Problem	Solution
The temperature visualisation LED (3) does not light up.	Check the USB connection cable (7) is fully inserted in the USB-A connection port on the power bank.
Heating is disabled before the automatic switch-off function is activated (90 minutes).	Check the power bank is charged.
The product has been improperly used, or is damaged.	Take the heating pad to an authorised technical service centre to be checked.
The USB connection cable (7) or another internal wire is damaged.	Take the heating pad to an authorised technical service centre to be checked.
There are holes or tears in the fabric, or definite folds that won't open out.	Take the heating pad to an authorised technical service centre to be checked.

DISPOSAL

The packaging of the appliance is made of recyclable materials. Dispose of the packaging in accordance with the environmental protection regulations.

Pursuant to Directive 2012/19/EU on waste electrical and electronic equipment (WEEE), the crossed bin symbol on the appliance or its packaging indicates that the appliance must be disposed of separately from other waste at the end of its useful life.

The user must therefore take the appliance to an appropriate recycling centre for electrical and electronic equipment. Alternatively, instead of handling the disposal, the user can hand over the appliance to the retailer when a new, equivalent one is purchased. The user can hand over electronic products measuring less than 25 cm to retailers of electronic products having a sales area of at least 400 m² free of charge and with no obligation to buy. Proper waste sorting prior to sending the equipment no longer used for recycling, processing and disposal in a manner that respects the environment contributes towards avoiding negative impacts on the environment and on health and promotes the reuse and/or recycling of the equipment's material.

SUPPORT AND WARRANTY

The appliance is guaranteed for a period of two years from the date of delivery. The warranty period starts from the date indicated on the receipt/invoice (if clearly legible), unless the purchaser can demonstrate that it was delivered at a later date. In the case of a product defect that existed prior to the delivery date, its repair or replacement is guaranteed without any cost, unless one of

the two remedies is disproportionate with respect to the other. The purchaser is required to report the lack of conformity to an authorised Service Centre within two months of the discovery of the defect.

The warranty does not cover all the parts that are defective as a result of the following:


- a. damage from transport or accidental drops,
- b. incorrect installation or inadequate electrical system,
- c. repairs or alterations performed by unauthorised personnel,
- d. insufficient or improper maintenance and cleaning,
- e. products and/or product parts subject to wear and/or consumables, including the reduced operating time of the batteries (if installed) due to usage or age,
- f. non-observance of the operating instructions for the appliance, negligent or careless usage,
- g. cable twisted and/or bent excessively,
- h. cable pulled or torn,
- i. heating panel bent during use, crumpled and/or torn,
- j. incorrect washing (e.g. centrifuge, wringing, with control unit enabled, etc.),
- k. obvious burn marks from external sources (cigarette etc.),
- l. bent back on itself.

The above list is provided as a mere example, and is not complete, as this warranty does not apply to any of the circumstances that cannot be connected to appliance manufacturing defects. The warranty is also excluded in the case of improper use of the appliance and in case of professional use. The manufacturer declines all liability for any damage that can be caused directly or indirectly to people, property or pets as a result of the failure to observe all the instructions indicated in the "Instructions and warnings booklet"

regarding appliance installation, use and maintenance. This does not apply to the seller's contractual warranty actions.

Service mode

The appliance may only be repaired by an authorised service centre. If under warranty, the defective appliance must be sent to the Service Centre together with proof of purchase that indicates the date of purchase or delivery.

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Tenacta Group S.p.A. Via Piemonte 5 IT-24052 Azzano S. P.